



# Deft Reactor™

The Digital Enabled  
Formulation Toolset Reactor

Your Instruments  
Your Workflow  
Made **Digital**


Deft  
Reactor



## Connected Digital Workflow

Working with partners GPE Scientific and major instrument manufacturer Julabo, the new **Deft Reactor** system brings device integration, digital control and reporting features of deft to lab reactors.

Common lab instruments such as the Julabo PRESTO thermo-circulator, overhead stirrers and pH meters are integrated with jacketed and unjacketed reactors such as Chemglass, to develop a highly configurable reactor that can be used either as standalone or alongside an existing Deft system.



## Your Workflow

The connected lab instruments via the **Deft Reactor**, can be fully controlled via the user-friendly, intuitive, dedicated App. Data can be input manually or automatically from a connected LIMS or ERP system putting you in charge of your experimental design and workflow.

You can operate and view the progress of your experiment remotely in real-time via the App.

Over-the-Air updates ensure your system is always state of the art and easily adaptable to the latest techniques.

